

# 4410/4412 DDS Urukul

### Features

- 4-channel 1GS/s DDS
- Output frequency from <1 to >400 MHz
- Sub-Hz frequency resolution
- Controlled phase steps
- Accurate output amplitude control

# Applications

- Dynamic low-noise RF source
- Driving RF electrodes in ion traps
- Driving acousto-optic modulators
- Form a laser intensity servo with 5108 Sampler

# **General Description**

The 4410/4412 DDS Urukul card is a 4hp EEM module, part of the ARTIQ/Sinara family. It adds frequency generation capabilities to carrier cards such as 1124 Kasli and 1125 Kasli-SoC. It can also be combined with 5018 ADC Sampler to form the ARTIQ SU-Servo configuration.

It provides 4 channels of DDS (direct digital synthesis) at 1GS/s. Output frequencies from <1 to >400 MHz are supported. The nominal maximum output power of each channel is 10dBm. Each channel can be attenuated from 0 to -31.5 dB by a digital attenuator. RF switches (1ns temporal resolution) on each channel provide 70 dB isolation.

4410 DDS Urukul features AD9910 chips, while 4412 DDS Urukul features AD9912 chips. These offer slightly different features and specifications. See below for a comparison between the two.



Figure 1: Simplified Block Diagram



Figure 2: Simplified DDS Signal Path



Figure 3: Urukul card and front panel

### Source

4410/4412 DDS Urukul, like all the Sinara hardware family, is open-source hardware, and design files (schematics, PCB layouts, BOMs) can be found in detail at the repository https://github.com/sinara-hw/Urukul/.

### Model comparison

4410 DDS Urukul uses AD9910<sup>1</sup> chips, whereas 4412 DDS Urukul uses AD9912<sup>2</sup> chips. In general, 4412/AD9912 is capable of much higher frequency resolution, at the cost of more detailed control features provided by 4410/AD9910, such as phase synchronization, digital ramp modulation or DRG, and digital amplitude control. See individual DDS IC datasheets for feature details, especially of AD9910, or ARTIQ code section below for examples of more complex experiments only possible with AD9910. The SUServo configuration is only available for AD9910.

Parameter	AD9910	AD9912	Unit	Conditions
Resolution				
Frequency <sup>1,5</sup>	0.25		Hz	
		8	μHz	
Phase offset <sup>1,2</sup>	16	14	bits	
DAC full scale current <sup>1,2</sup>	8	10	bits	
Digital amplitude <sup>1</sup>	14		bits	
Power consumption <sup>5</sup>	7	6.5	W	4x 400 MHz, 10.5 dBm, 52°C
ARTIQ SUServo available	Yes	No		

 Table 1: Specification Differences

# **Electrical Specifications**

Specifications of parameters are based on the datasheets of the DDS IC (AD9910<sup>1</sup>, AD9912AD9912<sup>2</sup>), clock buffer IC (Si53312<sup>3</sup>), digital attenuator IC (HMC542BLP4E<sup>4</sup>), Sinara project information<sup>5</sup> and corresponding test results<sup>6</sup>.

Parameter	Min.	Тур.	Max.	Unit	Conditions
Clock input					
Input frequency <sup>1,2</sup>	10		1000	MHz	PLL disabled
	3.2		60	MHz	AD9910, PLL enabled, no clock division
	12.8		240	MHz	AD9910, PLL enabled, 4x clock division
	11		200	MHz	AD9912, PLL enabled, no clock division
	44		800	MHz	AD9912, PLL enabled, 4x clock division

 Table 2: Recommended Operating Conditions

<sup>1</sup>https://www.analog.com/media/en/technical-documentation/data-sheets/AD9910.pdf

<sup>2</sup>https://www.analog.com/media/en/technical-documentation/data-sheets/AD9912.pdf

<sup>3</sup>https://www.skyworksinc.com/-/media/SkyWorks/SL/documents/public/data-sheets/Si5331x\_datasheet.pdf

<sup>4</sup>https://www.analog.com/media/en/technical-documentation/data-sheets/hmc542b.pdf

<sup>5</sup>https://github.com/sinara-hw/Urukul/wiki#details-specification-and-typical-performance-data

<sup>6</sup>https://github.com/sinara-hw/sinara/issues/354#issuecomment-352859041

### Table 3: Recommended Operating Conditions, cont.

Parameter	Min.	Тур.	Max.	Unit	Conditions
Nominal input power <sup>3</sup>		10		dBm	

### Table 4: RF Output Specifications

Parameter	Min.	Тур.	Max.	Unit	Conditions
Low frequency power <sup>6</sup>			-20	dBm	100 kHz output
			10	dBm	1 MHz output
Frequency <sup>5</sup>	1		400	MHz	
Digital attenuation <sup>4</sup>	-31.5		0	dB	
Resolution					
Temporal (I/O Update) <sup>5</sup>		4		ns	
Digital attenuation <sup>4</sup>		0.5		dB	

#### **Table 5: Electrical Characteristics**

Parameter	Min.	Тур.	Max.	Unit	Conditions
Digital attenuator glitch duration <sup>6</sup>		100		ns	
RF switch <sup>6</sup>					
Rise to 90%		100		ns	
Isolation		70		dB	
Turn-on chirp			0.1	deg/s	Excluding the first $\mu$ s
Crosstalk <sup>6</sup>		-84		dB	Victim RF switch opened
		-110		dB	Victim RF switch closed
Cross-channel-intermodulation <sup>6</sup>		-90		dB	
Phase noise <sup>6</sup>		-85		dBc/Hz	0.1 Hz
		-95		dBc/Hz	1 Hz
		-107		dBc/Hz	10 Hz
		-116		dBc/Hz	100 Hz
		-126		dBc/Hz	1 kHz
		-133		dBc/Hz	10 kHz
		-135		dBc/Hz	100 kHz
		-128		dBc/Hz	1 MHz
		-149		dBc/Hz	10 MHz

Parameter	Min.	Тур.	Max.	Unit	Conditions
Second-order harmonics <sup>6</sup>		-40		dB	6 dBm output
		-34		dB	10.5 dBm output
Third-order harmonics <sup>6</sup>		-54		dB	6 dBm output
		-28		dB	10.5 dBm output

Table 6: Electrical Characteristics, cont.

The tabulated performance characteristics above were produced using the following setup unless otherwise noted:

- 100 MHz input clock into SMA, 10 dBm
- Input clock divided by 4
- PLL with x40 multiplier
- Output frequency at 80 MHz or 81 MHz

Harmonic content of the DDS signals from 4410 DDS Urukul is tabulated below<sup>7</sup>. An external 125 MHz clock signal was supplied.

Frequency (MHz)	Output power (dBm) of the n <sup>th</sup> -order harmonic											
riequency (willz)	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>	9 <sup>th</sup>			
0.1	-21.14	-59.03	-54.93	-93.07	-73.38	-94.07	-84.78	-91.77	-96.61			
0.5	4.51	-15.45	-11.61	-25.02	-24.35	-51.70	-35.14	-34.46	-37.85			
1	7.67	-16.80	-12.32	-18.27	-29.25	-30.87	-34.51	-39.28	-39.84			
10	10.67	-12.69	-13.94	-26.12	-27.76	-36.11	-55.32	-43.85	-42.65			
20	10.86	-24.90	-13.65	-22.87	-28.67	-47.68	-35.85	-35.45	-38.48			
50	10.74	-14.18	-15.01	-27.57	-29.01	-38.05	-51.52	-44.53	-42.71			
100	9.70	-33.59	-16.72	-34.36	-26.81	-40.14	-41.07	-43.88	-56.89			
200	8.97	-22.22	-16.23	-24.89	-30.49	-37.97	-37.79	-38.80	-40.14			
300	8.27	-19.17	-19.51	-29.80	-34.75	-38.90	-51.92	-53.38	-57.95			
400	7.68	-17.82	-21.60	-33.04	-37.80	-50.37	-57.45	-59.80	-64.68			
500	-1.80	-41.57	-51.71	-72.36	-89.35	-91.63	-93.15	-84.54	-107.57			

Table 7: Harmonic content with 0.0 dB digital attenuation

7https://github.com/sinara-hw/Urukul/issues/29

Frequency (MHz)			Output	power (dB	m) of the	n <sup>th</sup> -order h	narmonic		
	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>	9 <sup>th</sup>
0.1	-27.06	-81.35	-62.09	-97.37	-84.11	-103.78	-91.37	-100.48	-104.22
0.5	-3.2	-37.82	-52.21	-66.76	-77.86	-85.92	-86.37	-97.59	-120.76
1	-0.43	-34.47	-47.80	-75.28	-86.45	-101.91	-93.22	-96.14	-106.71
10	1.95	-31.04	-28.23	-51.76	-57.29	-76.26	-78.15	-83.85	-80.20
20	2.10	-33.05	-28.30	-54.50	-52.31	-72.39	-70.96	-82.98	-82.58
50	1.89	-33.24	-28.50	-52.67	-48.35	-74.77	-77.26	-79.33	-73.58
100	0.80	-38.51	-63.22	-61.73	-71.97	-97.45	-97.67	-107.40	-93.03
200	0.05	-38.25	-42.16	-63.01	-84.55	-82.66	-108.85	-116.62	-99.45
300	-0.51	-35.91	-48.83	-82.43	-100.53	-111.79	-118.62	-120.05	-97.72
400	-1.20	-38.37	-49.77	-89.45	-74.66	-108.12	-116.75	-114.08	-102.29
500	-11.20	-61.47	-77.59	-74.73	-100.23	-93.12	-99.83	-86.71	-112.63

### Table 8: Harmonic content with 10.0 dB digital attenuation

### Table 9: Harmonic content with 20.0 dB digital attenuation

Fraguancy (MHz)			Output	power (dB	m) of the	n <sup>th</sup> -order h	armonic		
riequency (will2)	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>	9 <sup>th</sup>
0.1	-31.06	-82.29	-68.34	-109.04	-92.48	-111.23	-99.94	-109.85	-112.36
0.5	-11.99	-56.69	-71.73	-95.76	-101.86	-114.37	-102.81	-106.94	-116.72
1	-9.94	-54.54	-56.49	-89.12	-105.94	-110.93	-102.79	-107.01	-117.29
10	-7.89	-50.19	-57.35	-91.36	-97.88	-107.95	-103.53	-96.04	-108.26
20	-7.79	-52.72	-58.03	-90.75	-99.82	-102.07	-101.55	-104.73	-103.31
50	-7.96	-52.36	-59.26	-84.44	-87.55	-86.88	-97.76	-92.61	-83.19
100	-9.04	-57.40	-61.76	-78.50	-91.80	-117.64	-107.40	-112.64	-102.07
200	-9.73	-57.39	-72.31	-72.66	-93.26	-95.95	-125.22	-122.35	-130.24
300	-10.27	-58.65	-74.60	-109.24	-107.74	-115.75	-125.36	-124.54	-98.86
400	-10.94	-59.62	-79.36	-98.48	-74.72	-111.95	-119.18	-114.63	-104.34
500	-21.00	-78.52	-99.07	-74.91	-99.55	-92.91	-103.02	-87.33	-114.87

Frequency (MHz)		Output power (dBm) of the n <sup>th</sup> -order harmonic											
Trequency (WITZ)	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>	9 <sup>th</sup>				
0.1	-37.89	-85.04	-77.41	-122.04	-114.29	-115.58	-110.65	-120.06	-123.70				
0.5	-22.38	-71.24	-89.84	-107.81	-108.76	-127.83	-114.12	-118.34	-127.07				
1	-21.01	-72.10	-90.08	-111.97	-111.30	-127.43	-114.38	-118.07	-128.06				
10	-19.22	-72.13	-90.74	-110.14	-105.28	-114.04	-113.51	-94.85	-116.15				
20	-19.28	-75.95	-94.72	-91.71	-107.55	-112.85	-112.24	-116.33	-114.02				
50	-19.27	-74.93	-92.21	-95.77	-101.06	-97.92	-108.30	-103.60	-93.96				
100	-20.27	-79.05	-87.48	-89.73	-104.00	-117.98	-112.12	-110.51	-105.80				
200	-21.19	-78.33	-106.81	-82.70	-92.31	-109.93	-133.86	-120.94	-102.95				
300	-21.58	-80.96	-112.44	-110.40	-108.11	-115.68	-122.51	-125.25	-99.63				
400	-22.44	-82.73	-105.55	-98.03	-74.84	-113.93	-119.41	-114.93	-104.55				
500	-31.73	-93.37	-99.74	-75.03	-99.27	-92.84	-104.14	-87.46	-116.22				

#### Table 10: Harmonic content with 31.5 dB digital attenuation

The RMS voltage of a 4410 DDS Urukul channel at different amplitude scale factors is measured below. The DDS channel is directly connected to an oscilloscope with a  $50\Omega$ ~termination. The reported values are obtained from the oscilloscope.



Figure 4: RMS voltage, 0dB attenuation



The ideal RMS voltage is described by the linear function  $V_{\text{rms,ideal}}(\text{ASF}) = \frac{V_{\text{rms}}(0.1)}{0.1} * \text{ASF}$ . The measured RMS voltage divided by the full scale ideal RMS voltage (i.e.  $V_{\text{rms,ideal}}(1)$ ) is shown below.



Figure 6: RMS voltage scaled by ideal voltage at ASF=1, 100 MHz



Figure 7: Phase noise of Urukul clocked by internal oscillator



Figure 8: Phase noise of 200 MHz DDS Output



Figure 9: Harmonic content of 200 MHz DDS Output<sup>8</sup>



Figure 10: Harmonic content of 80 MHz DDS Output (6 dBm)<sup>6</sup>



Figure 11: Harmonic content of 80 MHz DDS Output (10 dBm)<sup>6</sup>



Figure 12: RF switch turn on transient<sup>6</sup>



Figure 13: Nyquist rejection 400 MHz to 600 MHz<sup>6</sup>

<sup>8</sup>https://github.com/sinara-hw/Urukul/issues/64



Figure 14: Nyquist rejection 450 MHz to 550 MHz<sup>6</sup>



Figure 15: Attenuator step from 20 to 60 digital (16+4dB switch glitch)<sup>6</sup>



Figure 16: Attenuator step from 31 to 32 digital (major carry glitch)<sup>6</sup>

# Front panel LEDs

4410/4412 Urukul features a number of indicator LEDs in the front panel. Each of four channel SMA connectors is accompanied by a green LED, used to indicate that RF output is enabled, and a red LED, which activates to indicate a DDS synchronization/PLL issue. Note that when bypassing PLL in ARTIQ (see below) LED may stay on.

Two additional LEDs indicate power good (green) and overtemperature (red).

# **Configuring Operation Mode**

Mode of operation is specified by a DIP switch. The DIP switch can be found at the top right corner of the card. The following table summarizes the required setting for each mode. ✓ indicates ON, while X indicates OFF.

Mode	DIP Switch							
Widde	1	2	3	4				
Default	×	×	×	×				
SU-Servo	1	1	×	×				

#### Table 11: DIP switch configurations



Figure 17: Position of DIP switch

# **Urukul Single-/Double-EEM Modes**

4410/4412 DDS Urukul cards can operate with either a single or double EEM connections. When only EEM0 is connected, the card will act in single-EEM mode; when both EEM0 and EEM1 are connected, the card will act in double-EEM mode.

Double-EEM mode additionally provides 1 ns temporal resolution RF switches and multi-chip synchronisation. Without EEM1, the only way to access the switches is through the CPLD, using SPI. With EEM1, RF switches can be controlled as a TTL output through the LVDS transceiver. 1 ns temporal resolution can then be achieved using the ARTIQ RTIO system. Double-EEM mode is also recommended for the SUServo configuration.

# **ARTIQ System Description Entry**

ARTIQ/Sinara firmware/gateware is generated according to a JSON system description file, allowing gateware to be specific to and optimized for a certain system configuration.

4410/4412 Urukul should be entered in the peripherals list of the corresponding core device in the following format:

```
{
"type": "urukul",
"dds": "ad9910", // or "ad9912", as appropriate
"ports": [0, 1], // second port is optional
"synchronization": true, // or false, for AD9910 only
"clk_sel": 2, // select 0 to 2 for clock source
"pll_en": 0 // PLL bypass, for higher external frequencies
"refclk": 125e6, // for external clock signal
}
```

Replace 0 and 1 with the EEM port numbers used on the core device. Any ports can be used. For single-EEM mode, simply specify only one port. The synchronization field is boolean, false by default, and only applies to AD9910. In the clk\_sel field, 0 represents the internal 100 MHz oscillator, 1 represents SMA input, and 2 represents MMCX input. The pll\_en field may be specified 0 or 1 and is 1 by default.

Note that the SUServo configuration requires a different system description entry. See SUServo section below.

# Example ARTIQ Code

The sections below demonstrate simple usage scenarios of extensions on the ARTIQ control system. These extensions make use of the resources of the 4410/4412 DDS Urukul. They do not exhaustively demonstrate all the features of the ARTIQ system.

The full documentation for ARTIQ software and gateware, including guides for their use, is available at https://m-labs.hk/artiq/manual/. Please consult the manual for details and reference material of the functions and structures used here.

For details of AD9910 capabilities, operation modes, profiles, signals, etc., see also the corresponding datasheet, e.g. https://www.analog.com/media/en/technical-documentation/data-sheets/AD9910.pdf.

### 10 MHz sinusoidal wave

Generates a 10MHz sinusoid from RF0 with full scale amplitude, attenuated by 6 dB. Both the CPLD and the DDS channels should be initialized. By default, AD9910 single-tone profiles are programmed to profile 7.

```
@kernel
def run(self):
    self.core.reset()
    self.cpld.init()
    self.dds0.init()
    self.dds0.cfg_sw(True)
    self.dds0.set_att(6.*dB)
    self.dds0.set(10*MHz)
```

If the synchronization feature of AD9910 is enabled, RF signal across different channels of the same Urukul can be synchronized. For example, phase-coherent RF signal can be produced on both channel 0 and channel 1 after configuring an appropriate phase mode.

```
@kernel
def run(self):
    self.core.reset()
    self.cpld.init()
    self.dds0.init()
    self.dds0.cfg_sw(True)
    self.dds0.set_phase_mode(PHASE_MODE_TRACKING)
    self.dds0.set_att(6.*dB)
    self.dds1.init()
    self.dds1.cfg_sw(True)
    self.dds1.set_phase_mode(PHASE_MODE_TRACKING)
    self.dds1.set_att(6.*dB)
    self.dds1.set_phase_mode(PHASE_MODE_TRACKING)
    self.dds1.set(frequency=10*MHz, phase=0.0)
    self.dds1.set(frequency=10*MHz, phase=0.25)  # 0.25 turns phase offset
```

Note that the phase difference between the 2 channels might not be exactly 0.25 turns, but it is a constant. It can be negated by adjusting the phase parameter.

### Periodic RF pulse (AD9910 Only)

This example demonstrates that the RF signal can be modulated by amplitude using the RAM modulation feature of the AD9910. By default, RAM profiles are programmed to profile 0.

```
def prepare(self):
    self.amp = [0.0, 0.0, 0.0, 0.7, 0.0, 0.7, 0.7] # Reversed Order
    self.asf_ram = [0] * len(self.amp)
@kernel
def init_dds(self, dds):
    self.core.break_realtime()
    dds.init()
    dds.set_att(6.*dB)
    dds.cfg_sw(True)
@kernel
def configure_ram_mode(self, dds):
    self.core.break_realtime()
    dds.set_cfr1(ram_enable=0)
    self.cpld.io_update.pulse_mu(8)
    self.cpld.set_profile(0)
                                 # Enable the corresponding RAM profile
                                 # Profile 0 is the default
    dds.set profile ram(start=0, end=len(self.asf ram)-1,
        step=250, profile=0, mode=RAM_MODE_CONT_RAMPUP)
    self.cpld.io_update.pulse_mu(8)
    dds.amplitude_to_ram(self.amp, self.asf_ram)
    dds.write_ram(self.asf_ram)
    self.core.break realtime()
    dds.set(frequency=5*MHz, ram_destination=RAM_DEST_ASF)
    # Pass osk_enable=1 to set_cfr1() if it is not an amplitude RAM
    dds.set_cfr1(ram_enable=1, ram_destination=RAM_DEST_ASF)
    self.cpld.io_update.pulse_mu(8)
@kernel
def run(self):
    self.core.reset()
    self.core.break_realtime()
    self.cpld.init()
    self.init dds(self.dds0)
    self.configure_ram_mode(self.dds0)
```

The generated RF output of the above example consists of the following features in sequence:

- 1. A 5 MHz RF pulse for 2 microseconds.
- 2. No signal for 1 microseconds.
- 3. A 5 MHz RF pulse for 1 microseconds.
- 4. No signal for 3 microseconds.
- 5. Go back to item 1.

The expected waveform is plotted on the following figure. Note that phase of the RF pulses may drift gradually. Urukul was operated with  $50\Omega$  termination for this waveform.



### Simple amplitude ramp (AD9910 only)

An amplitude ramp of an RF signal can be generated by modifying the self.amp array in the previous example.

```
def prepare(self):
    # Reversed Order
    self.amp = [1.0, 0.9, 0.8, 0.7, 0.6, 0.5, 0.4, 0.3, 0.2, 0.1, 0.0]
    self.asf_ram = [0] * len(self.amp)
```

The generated RF output has an incrementing amplitude scale factor (ASF), increasing by 0.1 at every microsecond. Once the ASF reaches 1.0, it drops back to 0.0 at the next microsecond. The expected waveform over 1 cycle is plotted on the following figure. Note that phase of the RF pulses may drift gradually. Urukul was operated with  $50\Omega$  termination for this waveform.



### RAM synchronization (AD9910 only)

Multiple RAM channels can also be synchronized. Similar to the 10 MHz single-tone RF signals, specify phase when calling dds.set() in configure\_ram\_mode. For example, set phase to 0 for the channels (phase=0.0):

```
dds.set(frequency=5*MHz, phase=0.0, ram_destination=RAM_DEST_ASF)
```

Then, replace the run() function with the following:

```
@kernel
def run(self):
    self.core.reset()
    self.core.break_realtime()
    self.cpld.init()
    self.init_dds(self.dds0)
    self.init_dds(self.dds1)
    self.dds0.set_phase_mode(PHASE_MODE_TRACKING)
    self.dds1.set_phase_mode(PHASE_MODE_TRACKING)
    self.configure_ram_mode(self.dds0)
    self.configure_ram_mode(self.dds1)
```

Two phase-coherent RF signals with the same waveform as the previous figure (from either RAM examples) should be generated.

### **ARTIQ SU-Servo**

ARTIQ also allows for the joint configuration of one or two 4410 DDS Urukul cards and a 5108 Sampler as an integrated servo for laser intensity stabilization and similar purposes. SU-Servo also supports other additional features, such as preconfigured profiles per channel and automatic integrator hold. Urukuls must be 4410 AD9910 variants (not 4412 AD9912) and set to SU-Servo mode by DIP switch.

### SU-Servo System Description Entry

SU-Servo should be entered in the peripherals list of the corresponding core device in a single entry in the following format:

```
{
    "type": "suservo",
    "sampler_ports": [0, 1],
    "urukul0_ports": [2, 3],
    "urukul1_ports": [4, 5], // optional
    "clk_sel": 2 // select 0 to 2
}
```

Enter the actual EEM port numbers used on the core device. Any ports can be used. If using only one Urukul (and half of the available Sampler channels), the urukull\_ports field may be left out entirely.

For the clk\_sel field, 0 represents the internal 100 MHz oscillator, 1 represents SMA input, and 2 represents MMCX input.

### Example SU-Servo Code

In SU-Servo configuration, amplitude of the Urukul DDS output can be controlled with the Sampler ADC input through PI control, characterised by the following transfer function:

$$H(s) = k_p + \frac{k_i}{s + \frac{k_i}{g}}$$

In the following example, the DDS amplitude is set proportionally to the ADC input from Sampler. We initialize SU-Servo and all channels first. Note that the programmable gain of the Sampler is  $10^0 = 1$  and the input range is [-10V, 10V].

```
@kernel
def run(self):
    self.core.reset()
    self.suservo.init()
    self.suservo.set_pgia_mu(0, 0) # unity gain
    self.suservo.cplds[0].set_att(0, 15.)
    self.suschannel0.set_y(profile=0, y=0.) # Clear integrator
```

Next, we set up the PI control as an IIR filter. It has -1 proportional gain  $k_p$  and no integrator gain  $k_i$ .

Then, configure the DDS frequency to 10 MHz with 3V input offset. When input voltage  $\geq$  offset voltage, the DDS output amplitude is 0.

```
self.suschannel0.set_dds(
    profile=0,
    offset=-.3, # 3 V with above PGIA settings
    frequency=10*MHz,
    phase=0.)
```

SU-Servo encodes the ADC voltage on a linear scale [-1, 1]. Therefore, 3V is converted to 0.3. Note that the ASF of all DDS channels ss capped at 1.0 and the amplitude clips when ADC input  $\leq -7V$  with the above IIR filter.

Finally, enable the SU-Servo channel with the IIR filter programmed beforehand:

```
self.suschannel0.set(en_out=1, en_iir=1, profile=0)
self.suservo.set_config(enable=1)
```

A 10 MHz DDS signal is generated from the example above, with amplitude controllable by ADC. The RMS voltage of the DDS channel against the ADC voltage is plotted. The DDS channel is terminated with  $50\Omega$ .



DDS signal should be attenuated. High output power affects the linearity due to the 1 dB compression point of the amplifier at 13 dBm output power. 15 dB attenuation at the digital attenuator was applied in this example.

# **Ordering Information**

To order, please visit https://m-labs.hk and choose 4410/4412 DDS Urukul in the ARTIQ/Sinara hardware selection tool. Cards can be ordered as part of a fully-featured ARTIQ/Sinara crate or standalone through the 'Spare cards' option. Otherwise, orders can also be made by writing directly to mailto:sales@m-labs.hk.

Information furnished by M-Labs Limited is provided in good faith in the hope that it will be useful. However, no responsibility is assumed by M-Labs Limited for its use. Specifications may be subject to change without notice.